



Insulated Shippers

Ship nonhazardous materials

- Add as many ice packs as you need—choose from a variety of sizes to meet your refrigeration needs

Purchase the size you need for efficient use of space and refrigerant. Compact shippers have a cylindrical interior to hold tubes, vials, and other small items. Medium- and large-capacity models have a rectangular interior and include cartons. Mailer 03726-53 has thick walls—save money on refrigerant. **Note:** For shipping nonhazardous materials only.

Catalog number	Inside Dimensions (L x W x H)		Outside Dimensions (L x W x H)		Wall thickness	Qty/ cs	Price/ ea
	in.	cm	in.	cm			
Compact insulated shippers							
DY-03730-20	3¼ H x 3 dia	8.3 H x 7.6 dia	6 x 6 x 6	15.2 x 15.2 x 15.2	1½"	12	
DY-03730-30	8 H x 4⅞ dia	20.3 H x 12.4 dia	8 x 8 x 11	20.3 x 20.3 x 27.9	1½"	6	
Medium-capacity shippers (L x W x H)							
DY-03726-11†	8 x 6 x 4	20.3 x 15.2 x 10.2	11 x 9 x 7	27.9 x 22.9 x 17.8	1½"	12	
DY-03726-21†	8 x 6 x 9	20.3 x 15.2 x 22.9	11 x 9 x 12	27.9 x 22.9 x 30.5	1½"	12	
DY-03726-31†	8 x 6 x 12	20.3 x 15.2 x 30.5	11 x 9 x 15	27.9 x 22.9 x 38.1	1½"	10	
DY-03726-53†	6½ x 5¾ x 8	16.5 x 14.6 x 20.3	11 x 10¼ x 12½	27.9 x 26.0 x 31.8	2¼"	8	
Large-capacity shippers (L x W x H)							
DY-03732-00‡	11¼ x 8½ x 11½	29.8 x 21.6 x 29.2	13¾ x 10½ x 13½	34.9 x 26.7 x 34.3	1"	8	
DY-03740-00‡	10¾ x 7½ x 10½	27.3 x 19.1 x 26.7	13¾ x 10½ x 13½	34.9 x 26.7 x 34.3	1½"	8	
DY-03736-01†	11½ x 9¾ x 6⅞	29.2 x 24.8 x 17.5	14⅞ x 13¼ x 10	37.8 x 33.7 x 25.4	1½"	8	
DY-03738-01†	11½ x 9¾ x 9¼	29.2 x 24.8 x 23.5	15½ x 13¾ x 12½	38.4 x 34.6 x 32.1	1½"	6	
DY-03739-00†	18¾ x 11⅞ x 12¼	47.6 x 30.2 x 31.1	22½ x 16¼ x 16½	56.2 x 41.3 x 41.0	1½"	2	

†Supplied with a 200-lb test corrugated carton. ‡Supplied with a 275-lb test fiberboard outer case.

Insulated Containers

A Tabletop Ice Storage Chest

- Holds up to 50 lb of dry ice

Tabletop chest has dry ice sublimation rate less than 2.5 lb per day at full 50 lb capacity. Chest is constructed from rugged polyethylene; molded bottom won't scratch tabletop. Inside dimensions are 11 1/2"L x 9 3/4"W x 9 1/4"H (29.2 x 24.8 x 23.5 cm). NSF certified.

B Insulated Field Carrier

- Handle for easy carrying—1 cu ft capacity

Insulated field carrier has efficient, 2 1/4"-thick polyurethane foam insulation core and rugged, double-wall polyethylene construction. Includes carrying handle, hinged lid, and adjustable latches. Inside dimensions are 16 5/8"L x 10 1/4"W x 10 3/4"H (42.2 x 26.0 x 27.3 cm). NSF certified.

C Large-Capacity Dry Ice Storage Chests

Chests contain extra-thick foam plastic insulation core. It takes six days for 90 lb of dry ice to sublimate to 5 lb in model 03721-02; 10 days for 155 lb in model 03721-00; and 15 days for 215 lb in model 03720-00. Models 03721-00 and 03721-02 have handles for easy carrying. NSF certified. Order dollies at right for improved mobility.



Key letter	Catalog number	Capacity		Inside dimensions (L x W x H)		Wall in. (cm)	Price
		cu ft	lb	in.	cm		
A	SN-03719-00	—	—	11 1/2 x 9 3/4 x 9 1/4	29.2 x 24.8 x 23.5	—	
B	SN-03720-20	—	—	16 5/8 x 10 1/4 x 10 3/4	42.2 x 26.0 x 27.3	—	
C	SN-03721-02	1.6	90	15 1/4 x 14 x 15 1/4	38.7 x 35.6 x 38.7	2 (5.1)	
	SN-03721-00	2.5	155	17 3/8 x 16 1/8 x 16 1/2	44.1 x 42.8 x 41.9	2 1/4 (6.4)	
	SN-03720-00	5	215	21 x 13 1/2 x 24	53.3 x 34.3 x 61.0	2 (5.1)	

Dollies for Large-Capacity Dry Ice Storage Chests

Move your large-capacity storage chests easily on these aluminum dollies, even when filled to capacity with dry ice. Rim on edge of platform keeps containers on dolly during acceleration or bump stops. Gray hammertone finish on heavy-gauge aluminum. Large 4" (10.2 cm) swivel casters roll easily. NSF certified.

SN-03721-10 Dolly for 03721-00 and -02 containers
SN-03720-10 Dolly for 03720-00 container